



Mata Kuliah: **Digital Signal Processing (IF3024)**

Tugas: **Final Project**

Daftar Anggota:

1. **Rustian Afencius Marbun (122140155)**
 2. **Shintya Ayu Wardani (122140138)**
-

Real Time Detection of rPPG and Respiratory Signals

1 Pendahuluan

1.1 Latar Belakang

Final project pada mata kuliah Pengolahan Sinyal Digital (IF3024) bertujuan untuk membangun sistem deteksi sinyal rPPG dan sinyal respirasi secara real-time menggunakan webcam. Sistem ini bekerja tanpa kontak fisik dengan memproses video langsung dari kamera untuk menampilkan estimasi detak jantung dan laju napas dalam antarmuka grafis.

Untuk sinyal rPPG, digunakan algoritma Plane Orthogonal-to-Skin (POS) yang menganalisis perubahan warna pada area wajah (dahi) untuk mengekstraksi sinyal denyut nadi. Metode POS dipilih karena terbukti baik pada pencahayaan dan gerakan ringan [1]. Sedangkan untuk sinyal respirasi, sistem mendeteksi gerakan bahu dari landmark pose MediaPipe. Perubahan posisi ini di filter menggunakan band-pass untuk menangkap frekuensi pernapasan (0.1–0.5 Hz)[2].

1.2 Tujuan

Tujuan dibuatnya sistem ini adalah sebagai berikut:

- Mengembangkan sistem pemantauan detak jantung dan laju napas secara real-time tanpa kontak menggunakan webcam.
- Mengimplementasikan algoritma POS untuk ekstraksi sinyal rPPG dan MediaPipe Pose + Butterworth filter untuk sinyal respirasi.
- Menyajikan hasil estimasi dalam bentuk grafik dan teks melalui antarmuka GUI yang interaktif.

1.3 Komponen Pendukung

Dalam penyelesaian final project ini, di perlukan alat dan bahan yang dapat mendukung dalam pembuatan program deteksi sinyal rPPG dan sinyal respirasi secara real time. Berikut komponen pendukung yang digunakan:

1.3.1 Bahasa Pemrograman

Dalam pengembangan program ini, bahasa pemrograman Python dipilih karena menyediakan library yang banyak digunakan dalam proyek berbasis computer vision.

1.3.2 Library

Berikut beberapa library yang digunakan dalam program ini :

1. **logging**: Mencatat log informasi, peringatan, dan kesalahan selama eksekusi aplikasi. Berguna untuk debugging dan pemantauan alur program.
2. **sys**: Mengakses parameter dan fungsi sistem. Digunakan untuk keluar dari aplikasi saat terjadi kesalahan kritis seperti kegagalan impor modul.
3. **tkinter**: Library GUI standar Python. Digunakan untuk membangun antarmuka pengguna seperti jendela utama, tombol, dan grafik.
4. **cv2 (OpenCV)**: Library visi komputer untuk pengambilan dan pemrosesan frame video dari kamera secara real-time, termasuk anotasi dan konversi warna.
5. **numpy**: Library komputasi numerik. Digunakan untuk manipulasi array, perhitungan statistik, dan pengolahan sinyal.
6. **matplotlib.pyplot**: Modul visualisasi data 2D. Digunakan untuk menampilkan grafik sinyal pernapasan dan detak jantung secara real-time.
7. **threading**: Mengelola eksekusi paralel (multithreading). Digunakan agar proses pengambilan frame tidak mengganggu GUI.
8. **queue**: Antrean thread-safe. Mengalirkan data frame antar-thread secara aman.
9. **time**: Mengatur delay, mengukur durasi, dan mencatat waktu pengambilan data.
10. **PIL (Pillow)**: Mengubah frame dari format OpenCV menjadi format gambar yang bisa ditampilkan di **tkinter**.
11. **collections.deque**: Buffer efisien untuk menyimpan dan membatasi jumlah sampel sinyal secara otomatis.
12. **mediapipe**: Framework real-time dari Google untuk deteksi pose dan wajah. Digunakan untuk mengekstrak landmark tubuh dan wajah.
13. **scipy**: Library ilmiah Python. Digunakan untuk pemrosesan sinyal lanjutan seperti filtering.
14. **pywt (PyWavelets)**: Digunakan untuk denoising sinyal menggunakan transformasi wavelet.
15. **dataclasses**: Modul Python untuk mendeklarasikan class data dengan sintaks yang ringkas dan bersih.
16. **typing**: Menyediakan petunjuk tipe (type hints) untuk meningkatkan keterbacaan dan pemeliharaan kode.
17. **asttokens, executing, stack-data, pure-eval**: Mendukung pelacakan eksekusi dan debugging stack trace.
18. **colorama, prompt-toolkit, pygments, wcwidth**: Untuk pewarnaan teks dan tampilan interaktif di terminal.
19. **comm, ipykernel, ipython, jupyter-client, jupyter-core, nest-asyncio, matplotlib-inline**: Untuk lingkungan Jupyter Notebook.
20. **debugpy**: Debugger Python yang digunakan oleh IDE seperti VSCode.

21. **decorator**: Untuk mempermudah pembuatan dekorator fungsi.
22. **exceptiongroup**: Untuk menangani banyak exception secara bersamaan.
23. **jedi, parso**: Untuk auto-completion dan parsing kode (biasanya digunakan di editor).
24. **packaging, platformdirs**: Untuk mengelola metadata paket dan direktori konfigurasi.
25. **psutil**: Untuk memantau resource sistem seperti CPU dan memori.
26. **pywin32**: Binding ke API Windows.
27. **pyzmq, tornado, traitlets**: Untuk komunikasi jaringan, digunakan dalam backend Jupyter.
28. **six**: Menjaga kompatibilitas lintas versi Python.
29. **typing-extensions**: Mendukung anotasi tipe baru pada Python versi lama.

2 Metode dan Analisis Matematis

2.1 Filter Respiration (Butterworth)

Untuk mengekstraksi sinyal respirasi, sistem ini menggunakan pergerakan vertikal bahu (koordinat Y dari landmark MediaPipe). Nilai ini direkam sebagai sinyal waktu, kemudian difilter menggunakan Butterworth band-pass filter orde 4 untuk mengisolasi frekuensi pernapasan manusia normal, yaitu sekitar 0.1–0.5 Hz (6–30 napas per menit)[2].

Fungsi Transfer

Filter Butterworth dipilih karena respon frekuensinya yang halus dan tanpa ripple pada passband maupun stopband[3]. Fungsi transfer analog dari Butterworth band-pass orde n secara umum adalah:

$$H(s) = \frac{1}{\sqrt{1 + \left(\frac{s}{\omega_c}\right)^{2n}}}$$

Untuk implementasi digital, digunakan fungsi desain dari `scipy.signal.butter()` yang mengubah frekuensi cutoff menjadi bentuk digital (diskret) berdasarkan frekuensi sampling f_s , dan diterapkan menggunakan zero-phase filtering dengan `scipy.signal.filtfilt()` untuk menghindari distorsi fase.

Parameter Filter:

- **Orde**: 4
- **Frek. sampling f_s** : 30 Hz (asumsi webcam 30 FPS)
- **Cutoff**:
 - Low = 0.1 Hz
 - High = 0.5 Hz

2.2 Algoritma rPPG (POS)

Untuk mendeteksi detak jantung, sistem ini menggunakan sinyal warna dari wajah (area dahi). Perubahan warna ini diakibatkan oleh fluktuasi volume darah akibat denyut jantung. Salah satu metode yang digunakan untuk mengekstraksi sinyal tersebut adalah Plane Orthogonal-to-Skin (POS)[1].

Metode POS memproyeksikan sinyal RGB dari video wajah ke arah yang ortogonal terhadap vektor warna kulit, sehingga menghasilkan sinyal dengan rasio sinyal terhadap noise (SNR) yang tinggi.

Langkah-langkah POS:

1. Ekstraksi sinyal RGB dari area wajah selama durasi tertentu.
2. Normalisasi dan pengurangan rata-rata per kanal (detrending warna):

$$\tilde{C}(t) = C(t) - \mu(C(t))$$

3. Proyeksi sinyal warna ke bidang ortogonal:


$$S(t) = P \cdot \tilde{C}(t)$$

4. Gabungkan dua sinyal hasil proyeksi dengan estimasi rasio energi untuk menghasilkan satu sinyal 1D rPPG yang dominan.
5. Estimasi detak jantung dari frekuensi dominan sinyal ini menggunakan FFT.

Kelebihan POS:

1. Robust terhadap variasi pencahayaan
2. Efektif dalam kondisi gerakan ringan
3. Lebih baik dibanding metode dasar seperti CHROM atau GREEN

3 Alur Pemrosesan



Gambar 1: Alur Sistem

4 Implementasi

4.1 Instalasi Library

4.1.1 Menggunakan Python venv

File konfigurasi yang ada di dalam direktori virtual environment Python yang dibuat. File ini bukanlah alat untuk menginstal library secara langsung.

```

1
2 python -m venv nama_env

```

Kode 1: Perintah create environment

4.1.2 Menggunakan Requirements.txt

file teks yang berisi daftar semua library Python dan versi spesifiknya yang dibutuhkan oleh suatu proyek.

```
1  
2 pip install -r requirements.txt
```

Kode 2: Perintah install requirements.txt

4.2 Folder Gui

Folder gui merupakan paket Python yang berisi komponen-komponen Graphical User Interface (GUI) untuk aplikasi pemantauan.

4.2.1 Modul pycache

- folder yang secara otomatis dibuat oleh Python ketika menjalankan script Python.
- Fungsi utamanya untuk menyimpan bytecode terkompilasi (.pyc file) dari modul-modul Python.
- Tujuannya untuk mempercepat waktu startup pada eksekusi selanjutnya, karena Python tidak perlu mengkompilasi ulang file sumber (.py file) setiap kali.

4.2.2 Modul init.py

Baris *UnifiedVitalSignsApp*, *PlotManager* mendefinisikan apa saja yang akan diekspor (dapat diakses) ketika paket gui diimpor menggunakan `from gui import *`.

```
1  
2 from gui.main_app import UnifiedVitalSignsApp  
3 from gui.plot_manager import PlotManager  
4  
5 __all__ = ['UnifiedVitalSignsApp', 'PlotManager']
```

Kode 3: init.py

4.2.3 Modul main-app.py

Modul ini berisi kelas utama aplikasi GUI, yaitu *UnifiedVitalSignsApp*.

1. Library

```
1  
2 import cv2  
3 import numpy as np  
4 import tkinter as tk  
5 from tkinter import ttk, messagebox  
6 import threading  
7 import queue  
8 import time  
9 from PIL import Image, ImageTk  
10 from typing import Optional, Tuple  
11  
12 from processors import RPPGProcessor, RespirationProcessor  
13 from gui.plot_manager import PlotManager  
14
```

Kode 4: library main-app.py

2. Kelas UnifiedVitalSignsApp mengoordinasikan GUI, pemrosesan video, dan visualisasi real-time untuk pemantauan detak jantung (menggunakan rPPG) dan laju pernapasan (menggunakan deteksi pose).

```

1
2 class UnifiedVitalSignsApp:
3     """
4     Unified GUI application for real-time vital signs monitoring.
5
6     This application provides a comprehensive interface for monitoring both
7     heart rate (using rPPG) and respiration rate (using pose detection)
8     from video input with real-time visualization and controls.
9     """
10
11     def __init__(self):
12         """Initialize the unified vital signs monitoring application."""
13         # Setup main window
14         self._setup_main_window()
15
16         # Initialize application state
17         self._initialize_application_state()
18
19         # Initialize signal processors
20         self._initialize_processors()
21
22         # Setup GUI components
23         self._setup_gui_components()
24
25         # Initialize threading components
26         self._initialize_threading()
27
28         # Bind window close event
29         self.root.protocol("WM_DELETE_WINDOW", self.on_closing)
30

```

Kode 5: Class UnifiedVitalSignsApp

3. Inisialisasi jendela utama aplikasi (setup-main-window).

```

1
2 def _setup_main_window(self) -> None:
3     """Setup the main application window with proper sizing and title."""
4     self.root = tk.Tk()
5     self.root.title("Heart Rate & Respiration Realtime Monitoring")
6     self.root.geometry("1400x900")
7     self.root.minsize(1200, 800) # Set minimum window size
8
9     # Configure main window style
10    style = ttk.Style()
11    style.theme_use('clam') # Use a modern theme
12
13    def _initialize_application_state(self) -> None:
14        """Initialize application state variables."""
15        self.is_running = False
16        self.cap: Optional[cv2.VideoCapture] = None
17        self.video_label: Optional[tk.Label] = None
18

```

Kode 6: Inisiasi Utama

4. Inisialisasi processor sinyal untuk rPPG dan pernapasan (initialize processors). Pengaturan komponen GUI seperti panel kontrol, tampilan video, dan area plot (setup-gui-components, setup-control-panel, setup-content-area, dll.).

```

1
2  def _initialize_processors(self) -> None:
3      """Initialize signal processing modules."""
4      # Initialize rPPG processor for heart rate detection
5      self.rppg_processor = RPPGProcessor(fps=30, window_length=1.6)
6
7      # Initialize respiration processor for breathing analysis
8      self.resp_processor = RespirationProcessor(fps=30)
9
10     def _setup_gui_components(self) -> None:
11         """Setup all GUI components and layout."""
12         # Create main container frame
13         main_frame = ttk.Frame(self.root)
14         main_frame.pack(fill=tk.BOTH, expand=True, padx=10, pady=10)
15
16         # Setup control panel
17         self._setup_control_panel(main_frame)
18
19         # Setup content area with video and plots
20         self._setup_content_area(main_frame)
21
22     def _setup_control_panel(self, parent: tk.Widget) -> None:
23         . . .
24
25     def _setup_content_area(self, parent: tk.Widget) -> None:
26         . . .
27
28     def _initialize_threading(self) -> None:
29         """Initialize threading components for video capture and processing."""
30         self.frame_queue = queue.Queue(maxsize=10)
31         self.capture_thread: Optional[threading.Thread] = None
32

```

Kode 7: Inisialisasi processor

5. Logika untuk memulai (start-monitoring) dan menghentikan (stop-monitoring) proses pemantauan, termasuk inisialisasi kamera dan pengelolaan thread.

```

1
2  def start_monitoring(self) -> None:
3      . . .
4  def stop_monitoring(self) -> None:
5      . . .
6      # Clear video display
7      cv2.destroyAllWindows()
8

```

Kode 8: Logika monitoring

6. Pengelolaan queue frame (capture-frames, process-frames) untuk memisahkan pengambilan video dari pemrosesan.

```

1
2  def _capture_frames(self) -> None:
3      . . .
4
5  def _process_frames(self) -> None:
6      . . .
7      # Schedule next frame processing
8      if self.is_running:
9          self.root.after(33, self._process_frames) # ~30 FPS

```


10

Kode 9: Pengelolaan frame

7. Memperbarui tampilan nilai vital sign (update-vital-sign-display) dan status aplikasi (update-status-display).

```

1
2     def _update_vital_signs_display(self, face_detected: bool, pose_detected: bool,
3                                     hr: float, rr: float) -> None:
4         . . .
5     def _update_status_display(self, face_detected: bool, pose_detected: bool) -> None:
6         . . .
7         self.plot_manager.update_respiration_data(
8             timestamp, raw_val, filtered_val
9         )
10

```

Kode 10: Tampilan vital sign

8. Mengelola tampilan frame video yang diproses (display-frame).

```

1
2     def _display_frame(self, frame: np.ndarray) -> None:
3         . . .
4
5     def _resize_frame_to_fit(self, frame: np.ndarray, width: int, height: int) -> np.ndarray:
6         . . .
7
8     return cv2.resize(frame, (new_width, new_height))
9

```

Kode 11: Tampilan frame video

9. Penanganan penutupan aplikasi (on-closing) untuk memastikan pelepasan sumber daya yang benar.

```

1
2     def on_closing(self) -> None:
3         """
4         Handle application closing event.
5
6         This method ensures proper cleanup of resources when the
7         application window is closed.
8         """
9         # Stop monitoring if running
10        self.stop_monitoring()
11
12        # Close any MediaPipe resources
13        if hasattr(self.resp_processor, 'landmarker') and self.resp_processor.landmarker:
14            self.resp_processor.landmarker.close()
15
16        # Destroy OpenCV windows
17        cv2.destroyAllWindows()
18
19        # Close tkinter application
20        self.root.quit()
21        self.root.destroy()
22
23    def run(self) -> None:
24        """
25        Start the unified vital signs monitoring application.

```

```

26
27     This method starts the main GUI event loop and runs the application
28     until the user closes the window.
29     """
30     try:
31         self.root.mainloop()
32     except KeyboardInterrupt:
33         print("Application interrupted by user")
34         self.on_closing()
35
36

```

Kode 12: Close aplikasi

4.2.4 Modul plot-manager.py

Modul yang mengelola semua fungsionalitas plotting Matplotlib untuk visualisasi data tanda-tanda vital secara real-time.

1. Library

```

1
2     import numpy as np
3     import matplotlib.pyplot as plt
4     from matplotlib.backends.backend_tkagg import FigureCanvasTkAgg
5     from matplotlib.animation import FuncAnimation
6     from collections import deque
7     import tkinter as tk
8     from typing import Tuple, List, Optional
9

```

Kode 13: Library plot-manager.py

2. Kelas PlotManager menangani pembuatan, pembaruan, dan animasi plot untuk rPPG dan sinyal pernapasan

```

1
2     class PlotManager:
3         """
4         Manages real-time plotting for rPPG and respiration signals.
5
6         This class handles the creation, updating, and animation of matplotlib
7         plots for visualizing heart rate and respiration monitoring data.
8         """
9
10        def __init__(self, rppg_parent: tk.Widget, resp_parent: tk.Widget):
11            """
12            Initialize the plot manager with parent widgets for embedding plots.
13
14            Args:
15                rppg_parent: Parent widget for rPPG plots
16                resp_parent: Parent widget for respiration plots
17            """
18            self.rppg_parent = rppg_parent
19            self.resp_parent = resp_parent
20
21            # Initialize plot data buffers
22            self._initialize_data_buffers()
23
24            # Setup plots
25            self._setup_rppg_plots()
26            self._setup_respiration_plots()

```

```

27
28         # Start animations
29         self._start_animations()
30

```

Kode 14: Library plot-manager.py

- Inisialisasi buffer data (initialize data buffers) menggunakan deque untuk menyimpan data sinyal terkini.
- Pengaturan plot rPPG (setup rppg plots) dengan tiga subplot (sinyal RGB, sinyal rPPG mentah, dan sinyal rPPG terfilter).
- Pengaturan plot pernapasan (setup respiration plots) dengan dua subplot (sinyal posisi bahu mentah dan sinyal pernapasan terfilter).
- Memulai animasi Matplotlib (start animations) menggunakan FuncAnimation untuk pembaruan real-time.
- Metode untuk memperbarui data plot rPPG (update rppg data) dan pernapasan (update respiration data).
- Fungsi pembaruan internal (update rppg plots, update resp plots) yang dipanggil oleh animasi untuk menggambar ulang grafik.
- Manajemen batas sumbu (update rppg axislimits, update resp axis limits) agar plot selalu menampilkan jendela waktu yang relevan dan penskalaan sumbu Y yang sesuai.
- Fungsionalitas untuk menghapus semua data plot (clear all data).

```

1
2     def _initialize_data_buffers(self) -> None:
3         . . .
4         # Respiration signal buffers
5         self.resp_time_data = deque(maxlen=buffer_size)
6         self.raw_resp_data = deque(maxlen=buffer_size)
7         self.filtered_resp_data = deque(maxlen=buffer_size)
8
9     def _setup_rppg_plots(self) -> None:
10        . . .
11
12    def _setup_filtered_rppg_subplot(self) -> None:
13        . . .
14
15    def _setup_filtered_respiration_subplot(self) -> None:
16        . . .
17        # Configure subplot appearance
18        ax.grid(True, alpha=0.3, linestyle='--')
19        ax.legend(loc='upper right', fontsize=9)
20        ax.set_facecolor('#f8f8f8')
21
22    def _start_animations(self) -> None:
23        . . .
24        # Animation for respiration plots (update every 100ms)
25        self.resp_ani = FuncAnimation(
26            self.resp_fig,
27            self._update_resp_plots,
28            interval=100,
29            blit=False,
30            cache_frame_data=False
31        )
32
33    def update_rppg_data(self, timestamp: float, r_val: float, g_val: float,

```

```

34         b_val: float, rppg_val: float, filtered_val: float) -> None:
35         . . .
36
37     def update_respiration_data(self, timestamp: float, raw_val: float,
38                               filtered_val: float) -> None:
39         . . .
40         self.resp_time_data.append(timestamp)
41         self.raw_resp_data.append(raw_val)
42         self.filtered_resp_data.append(filtered_val)
43
44     def _update_rppg_plots(self, frame) -> List:
45         . . .
46
47     def _update_rppg_axis_limits(self, time_array: np.ndarray) -> None:
48         . . .
49         except (IndexError, ValueError):
50             # Fallback to auto-scaling
51             ax.relim()
52             ax.autoscale_view()
53
54     def _update_resp_plots(self, frame) -> List:
55         . . .
56
57     def _update_resp_axis_limits(self, time_array: np.ndarray) -> None:
58         . . .
59         except (IndexError, ValueError):
60             ax.relim()
61             ax.autoscale_view()
62
63     def clear_all_data(self) -> None:
64         . . .
65         # Force plot updates
66         self._update_rppg_plots(None)
67         self._update_resp_plots(None)
68         self.rppg_canvas.draw_idle()
69         self.resp_canvas.draw_idle()
70
71     def save_plots(self, filename_prefix: str = "vital_signs") -> None:
72         """
73         Save current plots to image files.
74
75         Args:
76             filename_prefix: Prefix for saved filenames
77         """
78         try:
79             # Save rPPG plots
80             self.rppg_fig.savefig(f"{filename_prefix}_heart_rate.png",
81                                 dpi=300, bbox_inches='tight')
82
83             # Save respiration plots
84             self.resp_fig.savefig(f"{filename_prefix}_respiration.png",
85                                 dpi=300, bbox_inches='tight')
86
87             print(f"Plots saved as {filename_prefix}*.png")
88
89         except Exception as e:
90             print(f"Error saving plots: {e}")
91

```

Kode 15: inisialisasi

4.3 Folder Processors

Folder processors adalah paket Python yang berisi logika utama untuk pemrosesan sinyal vital, dibagi menjadi sub-paket untuk deteksi pernapasan dan deteksi detak jantung (rPPG).

4.3.1 Respiration

Modul ini mengimplementasikan deteksi laju pernapasan dengan menganalisis gerakan vertikal landmark bahu, yang mencerminkan ekspansi dan kontraksi dada selama siklus pernapasan. Ini menggunakan filter bandpass untuk mengisolasi frekuensi pernapasan yang relevan, biasanya dalam rentang 0.1 Hz hingga 0.5 Hz.

1. Modul Respiration-processor-modular.py

Modul ini merupakan versi yang di-refactor dari processor pernapasan, menggunakan komponen-komponen yang lebih kecil dan terfokus.

```

1
2     import numpy as np
3     import logging
4     from typing import Optional, Tuple
5     from utils.pose_detector import PoseDetectionHandler, PoseLandmarks
6     from utils.signal_buffer import SignalBufferManager, BufferConfig
7     from processors.respiration.respiratory_analyzer import RespiratorySignalAnalyzer,
8     FilterConfig, AnalysisResult
9     from utils.visualization_helper import VisualizationHelper, VisualizationConfig

```

Kode 16: Library Respiration-processor-modular.py

2. Class RespirationProcessorModular

- PoseDetectionHandler: Untuk mendeteksi landmark pose.
- SignalBufferManager: Untuk mengelola buffer sinyal.
- RespiratorySignalAnalyzer: Untuk menganalisis sinyal dan menghitung laju pernapasan.
- VisualizationHelper: Untuk menggambar visualisasi pada frame.

```

1
2     class RespirationProcessorModular:
3         """
4         Modular respiration processor using coordinated components.
5
6         This refactored version delegates specific responsibilities to focused
7         components, making the code more maintainable and testable.
8         """
9
10        def __init__(self, fps: int = 30, model_path: Optional[str] = None):
11            . . .
12            logging.info("Modular respiration processor initialized")
13
14        def _initialize_components(self, model_path: Optional[str]) -> None:
15            . . .
16            logging.info("All components initialized successfully")
17
18        def process_frame(self, frame: np.ndarray) -> Tuple[np.ndarray, bool, float]:
19            . . .
20
21        def _calculate_landmark_quality(self, landmarks: PoseLandmarks) -> float:
22            . . .

```

```

23         return np.clip(overall_quality, 0.0, 1.0)
24
25     def _analyze_respiratory_signal(self) -> None:
26         . . .
27         # Store filtered signal for visualization
28         if result.filtered_signal:
29             self.signal_buffer.add_filtered_sample(result.filtered_signal[-1])
30
31         logging.info(f"Respiration Rate: {self.current_rr:.1f} BPM (confidence: {
result.confidence:.3f})")
32         else:
33             logging.debug("Respiratory analysis failed or low confidence")
34
35     def update_filter_params(self, lowcut: float, highcut: float) -> None:
36         """
37         Update bandpass filter parameters.
38
39         Args:
40             lowcut: Low frequency cutoff in Hz
41             highcut: High frequency cutoff in Hz
42         """
43         self.signal_analyzer.update_filter_config(lowcut, highcut)
44         logging.info(f"Filter parameters updated: {lowcut:.2f} - {highcut:.2f} Hz")
45
46     def reset_signals(self) -> None:
47         """Reset all signal buffers and analysis state."""
48         self.signal_buffer.clear_buffers()
49         self.current_rr = 0
50         self.frame_idx = 0
51         self.last_analysis_result = None
52         logging.info("Respiration signals reset")
53
54     def get_signal_quality(self) -> float:
55         """
56         Get current signal quality metric.
57
58         Returns:
59             Signal quality score (0-1, higher is better)
60         """
61         return self.signal_buffer.get_overall_quality()
62
63     def get_buffer_statistics(self) -> dict:
64         . . .
65         return stats
66
67     def is_pose_detection_available(self) -> bool:
68         """
69         Check if pose detection is available.
70
71         Returns:
72             True if pose detection is working, False otherwise
73         """
74         return self.pose_detector.is_available()
75
76     def get_signal_data(self, num_samples: Optional[int] = None) -> dict:
77         """
78         Get current signal data for external analysis or visualization.
79
80         Args:
81             num_samples: Number of recent samples to retrieve
82
83         Returns:

```

```

84         Dictionary containing signal arrays
85         """
86         signal_data = self.signal_buffer.get_signal_data(num_samples)
87
88         return {
89             'raw_signal': signal_data.raw_signal,
90             'filtered_signal': signal_data.filtered_signal,
91             'timestamps': signal_data.timestamps,
92             'quality_scores': signal_data.quality_scores
93         }
94
95     def configure_visualization(self, **kwargs) -> None:
96         """
97         Configure visualization parameters.
98
99         Args:
100             **kwargs: Visualization configuration parameters
101         """
102         self.visualizer.update_config(**kwargs)
103         logging.info("Visualization configuration updated")
104
105     def validate_current_signal(self) -> Tuple[bool, str]:
106         """
107         Validate current signal data.
108
109         Returns:
110             Tuple of (is_valid, validation_message)
111         """
112         signal_data = self.signal_buffer.get_raw_signal()
113         return self.signal_analyzer.validate_signal(signal_data)
114

```

Kode 17: Class respiration modular

3. Modul Respiration-processor.py

Modul ini mengimplementasikan deteksi laju pernapasan dengan menganalisis gerakan vertikal landmark bahu.

```

1
2     import numpy as np
3     import mediapipe as mp
4     import cv2
5     import logging
6     from collections import deque
7     from typing import Optional, Tuple
8
9     from utils.signal_utils import (
10         apply_bandpass_filter,
11         find_dominant_frequency
12     )
13

```

Kode 18: Library Respiration-processor.py

4. Class RespirationProcessor

- Menginisialisasi sistem deteksi pose MediaPipe (setup pose detection), dengan opsi fallback ke deteksi pose dasar jika landmarker tingkat lanjut tidak tersedia.
- Mengelola buffer sinyal (initialize signal buffers) untuk menyimpan data pernapasan mentah dan yang sudah terfilter.

- Mengonfigurasi parameter filter bandpass (initialize filter parameters) untuk rentang frekuensi pernapasan (0.1 Hz - 0.5 Hz).
- Memproses setiap frame video untuk mengekstrak sinyal pernapasan, mendeteksi pose, dan menggambar visualisasi (process frame).
- Menghitung laju pernapasan dari gerakan bahu menggunakan analisis frekuensi (calculate respiration rate).
- Menyediakan fungsionalitas untuk memperbarui parameter filter (update filter params) dan mereset sinyal (reset signals).

```

1
2 class RespirationProcessor:
3     . . .
4
5     def __init__(self, fps: int = 30):
6         . . .
7         # Initialize respiration rate calculation
8         self.current_rr = 0
9         self.frame_idx = 0
10
11     def _setup_pose_detection(self) -> None:
12         . . .
13         self.pose_available = True
14         self.landmarker = None
15
16     def _initialize_signal_buffers(self) -> None:
17         """Initialize signal storage buffers for respiratory data."""
18         # Store up to 1 minute of data for analysis
19         self.max_buffer_size = self.fps * 60
20         self.raw_y_buffer = deque(maxlen=self.max_buffer_size)
21         self.filtered_y_buffer = deque(maxlen=self.max_buffer_size)
22         self.time_buffer = deque(maxlen=self.max_buffer_size)
23
24     def _initialize_filter_parameters(self) -> None:
25         """Initialize bandpass filter parameters for respiration frequencies."""
26         # Normal respiration rate: 12-20 breaths per minute (0.2-0.33 Hz)
27         # Allow wider range to capture individual variations
28         self.lowcut = 0.1 # Hz (6 breaths per minute)
29         self.highcut = 0.5 # Hz (30 breaths per minute)
30
31     def process_frame(self, frame: np.ndarray) -> Tuple[np.ndarray, bool, float]:
32         . . .
33
34     def _process_with_landmarker(self, frame: np.ndarray, processed_frame: np.ndarray,
35                                h: int, w: int) -> Tuple[bool, np.ndarray]:
36         . . .
37         pose_detected = False
38
39         # Convert to RGB and create MediaPipe image
40         rgb = cv2.cvtColor(frame, cv2.COLOR_BGR2RGB)
41         mp_image = mp.Image(image_format=mp.ImageFormat.SRGB, data=rgb)
42
43         # Calculate timestamp for video processing
44         timestamp_ms = int((self.frame_idx / self.fps) * 1000)
45         . . .
46         except Exception as e:
47             logging.error(f"Pose landmarker processing error: {e}")
48
49         return pose_detected, processed_frame
50
51     def _process_with_basic_pose(self, frame: np.ndarray, processed_frame: np.ndarray,

```



```

52         h: int, w: int) -> Tuple[bool, np.ndarray]:
53     . . .
54
55     def _draw_pose_landmarks(self, frame: np.ndarray, lx: int, ly: int,
56                             rx: int, ry: int) -> None:
57         . . .
58
59     def _store_respiratory_data(self, avg_y: float) -> None:
60         . . .
61
62     def _calculate_respiration_rate(self) -> None:
63         . . .
64         except Exception as e:
65             logging.error(f"Respiration calculation error: {e}")
66
67     def update_filter_params(self, lowcut: float, highcut: float) -> None:
68         """
69         Update bandpass filter parameters for respiration analysis.
70
71         Args:
72             lowcut: Low frequency cutoff in Hz
73             highcut: High frequency cutoff in Hz
74
75         Raises:
76             ValueError: If parameters are invalid
77         """
78         if lowcut >= highcut:
79             raise ValueError("Low cutoff must be less than high cutoff")
80
81         if lowcut < 0.05 or highcut > 1.0:
82             raise ValueError("Cutoff frequencies must be in range [0.05, 1.0] Hz")
83
84         self.lowcut = lowcut
85         self.highcut = highcut
86         logging.info(f"Updated filter parameters: {lowcut:.2f} - {highcut:.2f} Hz")
87
88     def reset_signals(self) -> None:
89         """Reset all signal buffers and respiration rate estimates."""
90         self.raw_y_buffer.clear()
91         self.filtered_y_buffer.clear()
92         self.time_buffer.clear()
93         self.current_rr = 0
94         self.frame_idx = 0
95         logging.info("Respiration signals reset")
96
97     def get_signal_quality(self) -> float:
98         """
99         Calculate and return current signal quality metric.
100
101         Returns:
102             Signal quality score (0-1, higher is better)
103         """
104         if len(self.raw_y_buffer) < 10:
105             return 0.0
106
107         # Calculate signal variability as quality indicator
108         recent_data = list(self.raw_y_buffer)[-30:] # Last 30 samples
109         signal_std = np.std(recent_data)
110         signal_range = np.max(recent_data) - np.min(recent_data)
111
112         # Normalize quality score (higher movement = better signal)
113         quality = min(1.0, signal_range / 50.0) # Assuming 50 pixels max movement

```

```

114
115         return quality
116
117

```

Kode 19: Class RespirationProcessor

5. Modul Respiration-analyzer.py

Modul ini fokus pada analisis sinyal pernapasan untuk perhitungan laju pernapasan.

```

1
2     import numpy as np
3     import logging
4     from typing import List, Optional, Tuple
5     from dataclasses import dataclass
6
7     from utils.signal_utils import (
8         apply_bandpass_filter,
9         find_dominant_frequency
10    )
11

```

Kode 20: Library Respiration-analyzer.py

6. Class RespiratorySignalAnalyzer

- Menerapkan filter bandpass untuk mengisolasi frekuensi pernapasan (filter signal).
- Menemukan frekuensi dominan dalam rentang pernapasan (find respiratory frequency).
- Menghitung kualitas sinyal (calculate signal quality) dan kepercayaan diri (calculate confidence) dalam hasil analisis.
- Menyediakan metode analyze signal yang mengembalikan objek AnalysisResult yang berisi laju pernapasan terhitung, frekuensi dominan, kualitas sinyal, dan kepercayaan diri.

```

1
2     @dataclass
3     class FilterConfig:
4         """
5         Configuration for respiratory signal filtering.
6
7         Attributes:
8             lowcut: Low frequency cutoff in Hz (minimum breathing rate)
9             highcut: High frequency cutoff in Hz (maximum breathing rate)
10            fs: Sampling frequency in Hz
11        """
12        lowcut: float = 0.1    # 6 breaths per minute
13        highcut: float = 0.5   # 30 breaths per minute
14        fs: float = 30.0       # 30 fps
15
16
17     @dataclass
18     class AnalysisResult:
19         """
20         Result of respiratory signal analysis.
21
22         Attributes:
23             respiration_rate: Calculated respiration rate in breaths per minute
24             dominant_frequency: Dominant frequency in Hz
25             signal_quality: Quality score of the analysis (0-1)
26             confidence: Confidence in the result (0-1)

```

```

27         filtered_signal: Filtered signal used for analysis
28         """
29         respiration_rate: float
30         dominant_frequency: float
31         signal_quality: float
32         confidence: float
33         filtered_signal: List[float]
34
35
36     class RespiratorySignalAnalyzer:
37         . . .
38
39         def __init__(self, filter_config: Optional[FilterConfig] = None):
40             . . .
41             logging.info(f"Respiratory analyzer initialized: {self.config.lowcut:.2f}-{self.
42             config.highcut:.2f} Hz")
43
44         def _validate_filter_config(self) -> None:
45             """Validate filter configuration parameters."""
46             if self.config.lowcut >= self.config.highcut:
47                 raise ValueError("Low cutoff must be less than high cutoff")
48
49             if self.config.lowcut < 0.05 or self.config.highcut > 1.0:
50                 raise ValueError("Cutoff frequencies must be in range [0.05, 1.0] Hz")
51
52             if self.config.fs <= 0:
53                 raise ValueError("Sampling frequency must be positive")
54
55         def analyze_signal(self, signal: List[float]) -> Optional[AnalysisResult]:
56             . . .
57             except Exception as e:
58                 logging.error(f"Respiratory analysis error: {e}")
59                 return None
60
61         def _filter_signal(self, signal: List[float]) -> List[float]:
62             . . .
63
64         def _find_respiratory_frequency(self, signal: List[float]) -> float:
65             . . .
66             except Exception as e:
67                 logging.error(f"Frequency analysis error: {e}")
68                 return 0.0
69
70         def _calculate_signal_quality(self, signal: List[float]) -> float:
71             . . .
72             # Calculate SNR
73             if noise_estimate > 0:
74                 snr = signal_power / noise_estimate
75                 # Normalize SNR to 0-1 range
76                 quality = min(1.0, snr / 10.0) # SNR of 10 = quality of 1.0
77             else:
78                 quality = 0.5 # Default if noise can't be estimated
79
80             return quality
81
82         def _calculate_confidence(self, signal: List[float], dominant_freq: float,
83             signal_quality: float) -> float:
84             . . . .
85             # Consider signal length (longer signals are more reliable)
86             length_factor = min(1.0, len(signal) / 100.0) # Full confidence at 100+ samples
87             confidence *= length_factor

```

```

88         # Ensure confidence is in valid range
89         return np.clip(confidence, 0.0, 1.0)
90
91     def update_filter_config(self, lowcut: float, highcut: float) -> None:
92         . . .
93
94     def get_filter_config(self) -> FilterConfig:
95         """
96         Get current filter configuration.
97
98         Returns:
99             Current FilterConfig object
100         """
101         return self.config
102
103     def validate_signal(self, signal: List[float]) -> Tuple[bool, str]:
104         . . .
105
106     def get_analysis_parameters(self) -> dict:
107         """
108         Get current analysis parameters.
109
110         Returns:
111             Dictionary containing analysis configuration
112         """
113         return {
114             'filter_lowcut': self.config.lowcut,
115             'filter_highcut': self.config.highcut,
116             'sampling_frequency': self.config.fs,
117             'min_analysis_length': self.min_analysis_length,
118             'confidence_threshold': self.confidence_threshold
119         }
120

```

Kode 21: Class RespiratorySignalAnalyzer

4.3.2 rPPG

1. Modul rppg-processor.py

Modul ini mengimplementasikan metode Plane-Orthogonal-to-Skin (POS) untuk mengekstraksi sinyal detak jantung dari video wajah.

```

1
2     import numpy as np
3     import mediapipe as mp
4     import cv2
5     import time
6     from collections import deque
7     from typing import Optional, Tuple, List
8
9     from utils.signal_utils import (
10         apply_bandpass_filter,
11         apply_savgol_filter,
12         wavelet_denoise,
13         normalize_signal,
14         detect_peaks_with_validation,
15         calculate_signal_quality,
16         smooth_signal_exponential
17     )

```

18

Kode 22: Library rppg-processor.py

2. Class RPPGProcessor

- Menginisialisasi deteksi wajah MediaPipe untuk mengidentifikasi ROI (Region of Interest) pada dahi (setup face detection).
- Mengekstrak ROI wajah (extract face roi).
- Mengelola buffer sinyal RGB dan sinyal rPPG yang diekstraksi (initialize signal buffers).
- Menerapkan metode POS (cpu pos) untuk memproyeksikan sinyal RGB ke bidang ortogonal terhadap warna kulit, meminimalkan artefak gerakan.
- Menerapkan berbagai filter dan pemrosesan sinyal (appl filter) seperti bandpass filtering, wavelet denoising, Savitzky-Golay smoothing, dan exponential smoothing.
- Menghitung detak jantung dari sinyal rPPG yang terfilter menggunakan deteksi puncak (calculate heart rate).
- Menyediakan metode untuk mereset sinyal (reset signals).

```

1
2  class RPPGProcessor:
3
4      def __init__(self, fps: int = 30, window_length: float = 1.6):
5          . . .
6          # Configure signal processing parameters
7          self._initialize_processing_parameters()
8
9      def _setup_face_detection(self) -> None:
10         """Initialize MediaPipe face detection with optimized settings."""
11         self.mp_face_detection = mp.solutions.face_detection
12         self.face_detection = self.mp_face_detection.FaceDetection(
13             model_selection=1, # Full-range model for better accuracy
14             min_detection_confidence=0.5
15         )
16
17     def _initialize_signal_buffers(self) -> None:
18         """Initialize all signal storage buffers."""
19         self.r_signal = deque(maxlen=self.max_buffer_size)
20         self.g_signal = deque(maxlen=self.max_buffer_size)
21         self.b_signal = deque(maxlen=self.max_buffer_size)
22         self.rppg_signal = deque(maxlen=self.max_buffer_size)
23         self.filtered_rppg = deque(maxlen=self.max_buffer_size)
24         self.timestamps = deque(maxlen=self.max_buffer_size)
25
26     def _initialize_hr_parameters(self) -> None:
27         """Initialize heart rate calculation parameters."""
28         self.current_hr = 0
29         self.hr_history = deque(maxlen=15) # Store recent HR estimates
30         self.hr_timestamps = deque(maxlen=15)
31         self.last_valid_hr_time = 0
32
33     def _initialize_processing_parameters(self) -> None:
34         """Initialize signal processing and filtering parameters."""
35         # Smoothing parameters
36         self.smooth_window_size = 9
37         self.wavelet_name = 'sym4'
38         self.wavelet_level = 3
39
40         # Peak detection parameters

```

```

41     self.min_peak_distance = int(self.fps * 0.5) # Minimum 0.5s between peaks
42     self.peak_prominence = 0.3
43     self.min_signal_quality = 1.2
44
45     # Timing control
46     self.force_recalc_interval = 1.0 # Force recalculation every second
47     self.last_calculation_time = 0
48     self.prev_filtered_value = 0
49     self.last_filter_time = 0
50
51     def cpu_pos(self, signal_array: np.ndarray) -> np.ndarray:
52         . . .
53         # Add to output with overlap-add
54         H[:, m:(n + 1)] = np.add(H[:, m:(n + 1)], Hnm)
55
56         return H[0, :]
57
58     def extract_face_roi(self, frame: np.ndarray) -> Tuple[Optional[np.ndarray], Optional
59 [Tuple[int, int, int, int]]]:
60         . . .
61         # Extract ROI
62         roi = frame[forehead_y:forehead_y+forehead_height,
63                     forehead_x:forehead_x+forehead_width]
64         bbox_coords = (forehead_x, forehead_y, forehead_width, forehead_height)
65
66         return roi, bbox_coords
67
68     def process_frame(self, frame: np.ndarray) -> Tuple[np.ndarray, bool, float]:
69         . . .
70         # Process signals if sufficient data available
71         if len(self.r_signal) >= self.window_size:
72             self._update_rppg_signal()
73             self._calculate_heart_rate()
74
75         return processed_frame, face_detected, self.current_hr
76
77     def _update_rppg_signal(self) -> None:
78         """Update rPPG signal using the POS method."""
79         if len(self.r_signal) < self.window_size:
80             return
81
82         # Prepare RGB array for POS processing
83         rgb_array = np.array([
84             list(self.r_signal),
85             list(self.g_signal),
86             list(self.b_signal)
87         ])
88
89         # Reshape for POS method: (batch_size=1, channels=3, time_samples)
90         rgb_array = rgb_array.reshape(1, 3, -1)
91
92         # Extract rPPG signal using POS method
93         rppg = self.cpu_pos(rgb_array)
94
95         # Store the latest rPPG value
96         if len(rppg) > 0:
97             self.rppg_signal.append(rppg[-1])
98
99         # Apply filtering if sufficient data available
100         if len(self.rppg_signal) >= self.window_size:
101             self._apply_filter()

```

```

102     def _apply_filter(self) -> None:
103         . . .
104         except Exception as e:
105             print(f"Filtering error: {e}")
106             # Fallback: use raw signal if filtering fails
107             if len(self.rppg_signal) > 0:
108                 self.filtered_rppg.append(self.rppg_signal[-1])
109
110     def _calculate_heart_rate(self) -> None:
111         . . .
112         except Exception as e:
113             print(f"Heart rate calculation error: {str(e)}")
114
115     def reset_signals(self) -> None:
116         """Reset all signal buffers and heart rate estimates."""
117         self.r_signal.clear()
118         self.g_signal.clear()
119         self.b_signal.clear()
120         self.rppg_signal.clear()
121         self.filtered_rppg.clear()
122         self.timestamps.clear()
123         self.hr_history.clear()
124         self.hr_timestamps.clear()
125         self.current_hr = 0
126         self.last_valid_hr_time = 0
127         self.prev_filtered_value = 0
128

```

Kode 23: Class RPPGProcessor

4.4 Folder Utils

Folder utils adalah paket Python yang berisi modul-modul utilitas atau pembantu yang menyediakan fungsionalitas umum yang dapat digunakan kembali di seluruh aplikasi. Ini membantu menjaga kode utama tetap bersih dan terfokus pada logika.

4.4.1 Modul pose-detector.py

Modul ini menyediakan handler khusus untuk pengaturan dan pemrosesan deteksi pose.

1. Library

```

1
2     import cv2
3     import numpy as np
4     import mediapipe as mp
5     import logging
6     from typing import Optional, Tuple, Union
7     from dataclasses import dataclass
8

```

Kode 24: Library pose-detector.py

2. Class PoseDetectionHandler

Class PoseDetectionHandler merangkum semua logika deteksi pose, mendukung MediaPipe Tasks yang lebih canggih dan estimasi pose dasar dengan fallback jika diperlukan.

- Menginisialisasi deteksi pose MediaPipe (setup pose detection), mencoba menggunakan landmarker tingkat lanjut terlebih dahulu dan beralih ke deteksi pose dasar jika terjadi kesalahan.

- Mendeteksi landmark pose dalam frame video (detect pose) dan mengembalikan objek PoseLandmarks yang berisi koordinat bahu dan confidence score.
- Memastikan apakah deteksi pose tersedia dan siap digunakan (is available).
- Objek PoseLandmarks digunakan sebagai dataclass untuk menyimpan data landmark pose yang terstruktur.

```

1
2 @dataclass
3 class PoseLandmarks:
4     """
5     Container for pose landmark data.
6
7     Attributes:
8         left_shoulder_x: X coordinate of left shoulder
9         left_shoulder_y: Y coordinate of left shoulder
10        right_shoulder_x: X coordinate of right shoulder
11        right_shoulder_y: Y coordinate of right shoulder
12        average_y: Average Y coordinate of both shoulders
13        confidence: Detection confidence score
14    """
15    left_shoulder_x: float
16    left_shoulder_y: float
17    right_shoulder_x: float
18    right_shoulder_y: float
19    average_y: float
20    confidence: float = 1.0
21
22    class PoseDetectionHandler:
23
24        def __init__(self, model_path: Optional[str] = None):
25            . . .
26
27        def _setup_pose_detection(self) -> None:
28            """
29            Initialize MediaPipe pose detection with fallback options.
30
31            Attempts to use the advanced pose landmarker if available,
32            falls back to basic pose estimation for compatibility.
33            """
34            try:
35                # Try to use advanced MediaPipe Tasks pose landmarker
36                self._setup_advanced_landmarker()
37                logging.info("Using advanced pose landmarker")
38
39            except Exception as e:
40                # Fallback to basic pose estimation
41                logging.info(f"Advanced pose landmarker not available ({e}), using basic pose
42            ")
43                self._setup_basic_pose()
44
45        def _setup_advanced_landmarker(self) -> None:
46            """Setup advanced MediaPipe Tasks pose landmarker."""
47            from mediapipe.tasks.python import BaseOptions
48            from mediapipe.tasks.python.vision import (
49                PoseLandmarker, PoseLandmarkerOptions, RunningMode
50            )
51
52            options = PoseLandmarkerOptions(
53                base_options=BaseOptions(model_asset_path=self.model_path),
54                running_mode=RunningMode.VIDEO,
55                num_poses=1

```



```

55         )
56         self.landmarker = PoseLandmarker.create_from_options(options)
57         self.pose_available = True
58
59     def _setup_basic_pose(self) -> None:
60         """Setup basic MediaPipe pose estimation."""
61         self.mp_pose = mp.solutions.pose
62         self.pose = self.mp_pose.Pose(
63             static_image_mode=False,
64             model_complexity=1,
65             smooth_landmarks=True,
66             min_detection_confidence=0.5,
67             min_tracking_confidence=0.5
68         )
69         self.pose_available = True
70         self.landmarker = None
71
72     def detect_pose(self, frame: np.ndarray, frame_idx: int, fps: float) -> Optional[
PoseLandmarks]:
73         . . .
74         if self.landmarker:
75             return self._detect_with_landmarker(frame, frame_idx, fps, h, w)
76         else:
77             return self._detect_with_basic_pose(frame, h, w)
78
79     def _detect_with_landmarker(self, frame: np.ndarray, frame_idx: int,
80                                fps: float, h: int, w: int) -> Optional[PoseLandmarks]:
81         . . .
82         except Exception as e:
83             logging.error(f"Pose landmarker processing error: {e}")
84
85         return None
86
87     def _detect_with_basic_pose(self, frame: np.ndarray, h: int, w: int) -> Optional[
PoseLandmarks]:
88         . . .
89         except Exception as e:
90             logging.error(f"Basic pose processing error: {e}")
91
92         return None
93
94     def is_available(self) -> bool:
95         """
96         Check if pose detection is available and ready.
97
98         Returns:
99             True if pose detection is available, False otherwise
100         """
101         return self.pose_available
102

```

Kode 25: Class PoseDetectionHandler

4.4.2 Modul signal-buffer.py

Modul ini menyediakan manajemen buffer sinyal yang efisien untuk pemrosesan sinyal pernapasan.

1. Library

```

1
2     import numpy as np
3     import logging

```

```

4 from collections import deque
5 from typing import List, Optional, Tuple
6 from dataclasses import dataclass, field
7

```

Kode 26: Library signal-buffer.py

2. Class SignalBufferManager

- Menginisialisasi buffer sinyal (deque) dengan ukuran maksimum yang ditentukan dalam BufferConfig.
- Menambahkan sampel baru ke buffer sinyal mentah, dengan timestamp dan quality score opsional (add sample).
- Menambahkan sampel sinyal yang sudah terfilter (add filtered sample).
- Mengambil data sinyal mentah, terfilter, timestamp, dan quality score (get raw signal, get filtered signal, get timestamps, get quality scores, get signal data).
- Memeriksa apakah buffer memiliki cukup data untuk analisis (has sufficient data).
- Mengosongkan semua buffer (clear buffers).
- Menghitung kualitas sinyal secara keseluruhan (get overall quality) dan statistik buffer (get buffer statistics).
- Dataclass BufferConfig dan SignalData digunakan untuk mengonfigurasi dan menyimpan data sinyal secara terstruktur.

```

1
2 @dataclass
3 class BufferConfig:
4     ...
5     max_size: int = 1800 # 1 minute at 30 fps
6     min_analysis_size: int = 30 # Minimum samples for analysis
7     fps: float = 30.0
8
9 @dataclass
10 class SignalData:
11     ...
12     raw_signal: List[float] = field(default_factory=list)
13     filtered_signal: List[float] = field(default_factory=list)
14     timestamps: List[float] = field(default_factory=list)
15     quality_scores: List[float] = field(default_factory=list)
16
17 class SignalBufferManager:
18     ...
19     logging.info(f"Signal buffer manager initialized with max size: {self.config.
20 max_size}")
21
22     def add_sample(self, value: float, quality_score: Optional[float] = None) -> None:
23         ...
24         logging.debug(f"Added sample: value={value:.2f}, timestamp={timestamp:.2f}s")
25
26     def add_filtered_sample(self, value: float) -> None:
27         """
28         Add a filtered signal sample.
29
30         Args:
31             value: Filtered signal value to add
32         """
33         self.filtered_buffer.append(value)

```

```

34     def get_raw_signal(self, num_samples: Optional[int] = None) -> List[float]:
35         . . .
36
37     def get_filtered_signal(self, num_samples: Optional[int] = None) -> List[float]:
38         . . .
39
40     def get_timestamps(self, num_samples: Optional[int] = None) -> List[float]:
41         . . .
42
43     def get_signal_data(self, num_samples: Optional[int] = None) -> SignalData:
44         . . .
45
46     def get_quality_scores(self, num_samples: Optional[int] = None) -> List[float]:
47         """
48         Get signal quality scores.
49
50         Args:
51             num_samples: Number of recent scores to retrieve (all if None)
52
53         Returns:
54             List of quality scores
55         """
56         if num_samples is None:
57             return list(self.quality_buffer)
58         else:
59             return list(self.quality_buffer)[-num_samples:] if len(self.quality_buffer)
60             >= num_samples else list(self.quality_buffer)
61
62     def has_sufficient_data(self) -> bool:
63         """
64         Check if buffer contains sufficient data for analysis.
65
66         Returns:
67             True if sufficient data is available, False otherwise
68         """
69         return len(self.raw_buffer) >= self.config.min_analysis_size
70
71     def get_buffer_size(self) -> int:
72         """
73         Get current buffer size.
74
75         Returns:
76             Number of samples currently in buffer
77         """
78         return len(self.raw_buffer)
79
80     def get_buffer_duration(self) -> float:
81         """
82         Get current buffer duration in seconds.
83
84         Returns:
85             Duration of buffered data in seconds
86         """
87         return len(self.raw_buffer) / self.config.fps
88
89     def clear_buffers(self) -> None:
90         """Clear all signal buffers and reset frame count."""
91         self.raw_buffer.clear()
92         self.filtered_buffer.clear()
93         self.time_buffer.clear()
94         self.quality_buffer.clear()
95         self.frame_count = 0

```

```

95         logging.info("Signal buffers cleared")
96
97     def get_overall_quality(self) -> float:
98         """
99         Calculate overall signal quality from recent samples.
100
101         Returns:
102             Overall quality score (0-1, higher is better)
103         """
104         if len(self.quality_buffer) == 0:
105             return 0.0
106
107         # Use recent samples for quality assessment
108         recent_quality = list(self.quality_buffer)[-30:] if len(self.quality_buffer) >=
30 else list(self.quality_buffer)
109         return np.mean(recent_quality)
110
111     def _calculate_sample_quality(self, value: float) -> float:
112         """
113         # Combine quality metrics
114         overall_quality = (range_quality + consistency_quality) / 2
115
116         return np.clip(overall_quality, 0.0, 1.0)
117
118     def get_buffer_statistics(self) -> dict:
119         """
120         return {
121             'size': len(raw_data),
122             'duration': self.get_buffer_duration(),
123             'mean': np.mean(raw_data),
124             'std': np.std(raw_data),
125             'range': np.max(raw_data) - np.min(raw_data),
126             'quality': self.get_overall_quality()
127         }
128
129

```

Kode 27: Class SignalBufferManager

4.4.3 Modul signal-utils.py

Modul ini berisi berbagai fungsi utilitas untuk pemrosesan sinyal.

```

1
2     import numpy as np
3     import scipy.signal as signal
4     import pywt
5     from typing import Tuple, List, Optional
6

```

Kode 28: Library signal-utils

2.
 - apply bandpass filter: Menerapkan filter bandpass ke sinyal untuk mengisolasi rentang frekuensi tertentu.
 - apply savgol filter: Menerapkan filter Savitzky Golay untuk penghalusan sinyal.
 - wavelet denoise: Melakukan denoising (pengurangan noise) menggunakan transformasi wavelet.
 - normalize signal: Normalisasi sinyal untuk konsistensi.
 - detect peaks with validation: Mendeteksi puncak dalam sinyal dengan validasi.

- calculate signal quality: Menghitung metrik kualitas sinyal.
- smooth signal exponential: Penghalusan sinyal menggunakan exponential smoothing.
- find dominant frequency: Menemukan frekuensi dominan dalam sinyal, kemungkinan menggunakan Fast Fourier Transform (FFT).

```

1
2 def apply_bandpass_filter(data: List[float], lowcut: float, highcut: float,
3                             fs: float, order: int = 4) -> np.ndarray:
4     ...
5     try:
6         b, a = signal.butter(order, [low, high], btype='band')
7         return signal.filtfilt(b, a, data)
8     except Exception as e:
9         raise ValueError(f"Filter design failed: {e}")
10
11
12 def apply_savgol_filter(data: np.ndarray, window_length: int = 15,
13                          polyorder: int = 2) -> np.ndarray:
14     ...
15     # Ensure window length is odd and valid
16     if window_length % 2 == 0:
17         window_length -= 1
18     window_length = max(3, min(window_length, len(data) - 2))
19
20     return signal.savgol_filter(data, window_length, polyorder)
21
22
23 def wavelet_denoise(data: np.ndarray, wavelet: str = 'sym4',
24                     levels: int = 3) -> np.ndarray:
25     ...
26     except Exception as e:
27         print(f"Wavelet denoising failed: {e}")
28         return data
29
30
31 def normalize_signal(signal_data: List[float]) -> np.ndarray:
32     ...
33     if std < 1e-9:
34         raise ValueError("Signal has zero variance - cannot normalize")
35
36     return (signal_array - mean) / std
37
38
39 def find_dominant_frequency(signal_data: np.ndarray, fs: float,
40                             freq_range: Optional[Tuple[float, float]] = None) -> float:
41     ...
42     # Find peak frequency
43     if len(fft_vals) == 0 or np.max(fft_vals) == 0:
44         return 0.0
45
46     peak_idx = np.argmax(fft_vals)
47     return freqs[peak_idx]
48
49
50 def detect_peaks_with_validation(signal_data: np.ndarray, fs: float,
51                                 min_distance_sec: float = 0.5,
52                                 prominence: float = 0.2,
53                                 min_width_sec: float = 0.08) -> Tuple[np.ndarray, List[
54 float]]:
55     ...
56     # Validate peak intervals for physiological plausibility
57     validated_rates = []

```

```

57     if len(peaks) >= 2:
58         for i in range(1, len(peaks)):
59             interval_samples = peaks[i] - peaks[i-1]
60             if interval_samples > 0:
61                 # Convert to rate (beats/breaths per minute)
62                 rate = 60.0 * fs / interval_samples
63                 # Accept rates in physiological range (40-180 BPM)
64                 if 40 <= rate <= 180:
65                     validated_rates.append(rate)
66
67     return peaks, validated_rates
68
69
70 def calculate_signal_quality(signal_data: np.ndarray) -> float:
71     . . .
72     # Combine amplitude and SNR for quality score
73     quality_score = amplitude_range * np.log10(snr + 1)
74     return max(0.0, quality_score)
75
76
77 def smooth_signal_exponential(current_value: float, previous_value: float,
78                               alpha: float = 0.3) -> float:
79     . . .
80     if not 0 < alpha < 1:
81         raise ValueError("Alpha must be between 0 and 1")
82
83     return alpha * current_value + (1 - alpha) * previous_value
84

```

Kode 29: Pemrosesan sinyal.

4.4.4 Modul visualization-helper.py

Modul ini bertanggung jawab untuk membantu visualisasi pada frame video yang diproses.

1. Library

```

1
2     import cv2
3     import numpy as np
4     from typing import Tuple, List, Optional
5     from dataclasses import dataclass
6
7     from utils.pose_detector import PoseLandmarks
8

```

Kode 30: Library visualization-helper.py

2. Class VisualizationHelper

- draw pose landmarks: Menggambar landmark pose pada frame.
- draw respiration rate: Menampilkan laju pernapasan pada frame.
- create signal overlay: Membuat overlay sinyal pada frame.
- Dataclass VisualizationConfig untuk mengonfigurasi parameter visualisasi.

```

1
2     @dataclass
3     class VisualizationConfig:
4         """
5         Configuration for visualization elements.

```

```

6
7     Attributes:
8         landmark_color: Color for pose landmarks (BGR format)
9         line_color: Color for connecting lines (BGR format)
10        text_color: Color for text labels (BGR format)
11        landmark_radius: Radius of landmark circles
12        line_thickness: Thickness of connecting lines
13        text_font: OpenCV font type
14        text_scale: Font scale factor
15        text_thickness: Text line thickness
16    """
17    landmark_color: Tuple[int, int, int] = (255, 0, 0) # Blue
18    line_color: Tuple[int, int, int] = (0, 255, 255) # Yellow
19    text_color: Tuple[int, int, int] = (255, 255, 255) # White
20    landmark_radius: int = 4
21    line_thickness: int = 2
22    text_font: int = cv2.FONT_HERSHEY_SIMPLEX
23    text_scale: float = 0.5
24    text_thickness: int = 1
25
26
27    class VisualizationHelper:
28
29        def __init__(self, config: Optional[VisualizationConfig] = None):
30            . . .
31            # Draw confidence indicator if available
32            if landmarks.confidence < 1.0:
33                self._draw_confidence_indicator(annotated_frame, landmarks.confidence)
34
35            return annotated_frame
36
37        def _draw_landmark(self, frame: np.ndarray, x: int, y: int, label: str) -> None:
38            """
39            Draw a single landmark point with label.
40
41            Args:
42                frame: Frame to draw on
43                x: X coordinate
44                y: Y coordinate
45                label: Text label for the landmark
46            """
47            # Draw landmark circle
48            cv2.circle(frame, (x, y), self.config.landmark_radius,
49                      self.config.landmark_color, -1)
50
51            # Draw label above the landmark
52            label_y = max(y - 10, 20) # Ensure label is visible
53            cv2.putText(frame, label, (x - 20, label_y),
54                      self.config.text_font, self.config.text_scale,
55                      self.config.text_color, self.config.text_thickness)
56
57        def _draw_roi_indicator(self, frame: np.ndarray, landmarks: PoseLandmarks) -> None:
58            """
59            Draw region of interest indicator.
60
61            Args:
62                frame: Frame to draw on
63                landmarks: Pose landmarks for ROI calculation
64            """
65            # Calculate ROI center
66            center_x = int((landmarks.left_shoulder_x + landmarks.right_shoulder_x) / 2)
67            center_y = int(landmarks.average_y)

```

```

68
69     # Draw ROI label
70     cv2.putText(frame, "Respiration ROI", (center_x - 50, center_y - 30),
71                 self.config.text_font, self.config.text_scale,
72                 self.config.text_color, self.config.text_thickness)
73
74     def _draw_confidence_indicator(self, frame: np.ndarray, confidence: float) -> None:
75         . . .
76         # Confidence text
77         cv2.putText(frame, f"Conf: {confidence:.2f}", (bar_x, bar_y - 5),
78                     self.config.text_font, self.config.text_scale,
79                     self.config.text_color, self.config.text_thickness)
80
81     def _get_confidence_color(self, confidence: float) -> Tuple[int, int, int]:
82         """
83         Get color based on confidence level.
84
85         Args:
86             confidence: Confidence score (0-1)
87
88         Returns:
89             BGR color tuple
90         """
91         if confidence > 0.8:
92             return (0, 255, 0)    # Green
93         elif confidence > 0.5:
94             return (0, 255, 255)  # Yellow
95         else:
96             return (0, 0, 255)    # Red
97
98     def draw_respiration_rate(self, frame: np.ndarray, respiration_rate: float,
99                             signal_quality: float = 1.0) -> np.ndarray:
100         . . .
101         # Draw quality indicator
102         self._draw_quality_indicator(annotated_frame, signal_quality, text_x, text_y -
25)
103
104         return annotated_frame
105
106     def _get_quality_color(self, quality: float) -> Tuple[int, int, int]:
107         """
108         Get color based on signal quality.
109
110         Args:
111             quality: Quality score (0-1)
112
113         Returns:
114             BGR color tuple
115         """
116         if quality > 0.7:
117             return (0, 255, 0)    # Green
118         elif quality > 0.4:
119             return (0, 255, 255)  # Yellow
120         else:
121             return (0, 128, 255)  # Orange
122
123     def _draw_quality_indicator(self, frame: np.ndarray, quality: float,
124                               x: int, y: int) -> None:
125         """
126         Draw signal quality indicator.
127
128         Args:

```



```

129         frame: Frame to draw on
130         quality: Quality score (0-1)
131         x: X position for indicator
132         y: Y position for indicator
133     """
134     quality_text = f"Quality: {quality:.2f}"
135     color = self._get_quality_color(quality)
136
137     cv2.putText(frame, quality_text, (x, y),
138                 self.config.text_font, self.config.text_scale,
139                 color, self.config.text_thickness)
140
141     def draw_breathing_indicator(self, frame: np.ndarray, breathing_phase: str,
142                                intensity: float = 0.5) -> np.ndarray:
143         ...
144         # Draw text
145         cv2.putText(annotated_frame, text, (indicator_x, indicator_y),
146                     self.config.text_font, self.config.text_scale * 1.2,
147                     color, self.config.text_thickness + 1)
148
149         return annotated_frame
150
151     def create_signal_overlay(self, frame: np.ndarray, signal_data: List[float],
152                              max_points: int = 100) -> np.ndarray:
153         ...
154         # Add overlay label
155         cv2.putText(annotated_frame, "Signal",
156                     (w - max_points * 2 - 15, overlay_y_start - 15),
157                     self.config.text_font, self.config.text_scale,
158                     self.config.text_color, self.config.text_thickness)
159
160         return annotated_frame
161
162     def update_config(self, **kwargs) -> None:
163         """
164         Update visualization configuration.
165
166         Args:
167             **kwargs: Configuration parameters to update
168         """
169         for key, value in kwargs.items():
170             if hasattr(self.config, key):
171                 setattr(self.config, key, value)
172             else:
173                 raise ValueError(f"Invalid configuration parameter: {key}")
174
175     def get_config(self) -> VisualizationConfig:
176         """
177         Get current visualization configuration.
178
179         Returns:
180             Current VisualizationConfig object
181         """
182         return self.config
183
184

```

Kode 31: Class

4.4.5 Modul Main.py

File main.py adalah entry point utama untuk aplikasi pemantauan tanda-tanda vital. Ini berfungsi sebagai orkestrator yang menginisialisasi sistem dan menjalankan aplikasi GUI.

1. Library

```

1
2     import logging
3     import sys
4     from tkinter import messagebox
5
6     # Import modular components
7     try:
8         from gui.main_app import UnifiedVitalSignsApp
9     except ImportError as e:
10         print(f"Error importing GUI components: {e}")
11         print("Please ensure all required packages are installed and modules are in the
12         correct location.")
13         sys.exit(1)

```

Kode 32: Library Main.py

2. Inisialisasi sistem

- Impor Modul: Mengimpor modul penting seperti logging, sys, tkinter.messagebox, dan UnifiedVitalSignsApp dari gui main app, serta menangani kegagalan impor.
- Fungsi setup_logging(): Mengonfigurasi pencatatan log ke konsol dan file, dengan level INFO dan format lengkap, serta mencatat pesan saat aplikasi dimulai.
- Fungsi check_dependencies(): Memverifikasi ketersediaan modul penting (cv2, numpy, scipy, dll). Jika ada yang hilang, tampilkan pesan kesalahan GUI dan log error.
- Fungsi main(): Menjalankan urutan inisialisasi aplikasi: logging, pengecekan dependency, peluncuran GUI (app.run()), dan penanganan kesalahan fatal secara terpusat.
- Blok ": Memastikan main() hanya dijalankan jika main.py dieksekusi langsung, bukan saat diimpor.

```

1
2     # Import modular components
3     try:
4         from gui.main_app import UnifiedVitalSignsApp
5     except ImportError as e:
6         print(f"Error importing GUI components: {e}")
7         print("Please ensure all required packages are installed and modules are in the
8         correct location.")
9         sys.exit(1)
10
11
12     def setup_logging():
13         """
14         Configure logging for application monitoring.
15
16         Sets up both console and file logging with appropriate formatting
17         to track application events and potential issues.
18         """
19         logging.basicConfig(
20             level=logging.INFO,
21             format="[%asctime)s] %(levelname)s - %(message)s",
22             handlers=[

```

```

22         logging.StreamHandler(),
23         logging.FileHandler('vital_signs_monitor.log')
24     ]
25 )
26
27 # Log application startup
28 logging.info("Vital Signs Monitor application starting...")
29
30
31 def check_dependencies():
32     """
33     Check if all required dependencies are available.
34
35     Returns:
36         bool: True if all dependencies are available, False otherwise
37     """
38     required_modules = [
39         'cv2', 'numpy', 'scipy', 'matplotlib', 'tkinter',
40         'mediapipe', 'PIL', 'pywt'
41     ]
42
43     missing_modules = []
44
45     for module in required_modules:
46         try:
47             __import__(module)
48         except ImportError:
49             missing_modules.append(module)
50
51     if missing_modules:
52         error_msg = f"Missing required modules: {' '.join(missing_modules)}"
53         logging.error(error_msg)
54         messagebox.showerror(
55             "Missing Dependencies",
56             f"{error_msg}\n\nPlease install the required packages using:\n\npip install -r requirements.txt"
57         )
58         return False
59
60     return True
61
62
63 def main():
64     """
65     Main entry point for the Vital Signs Monitor application.
66
67     This function:
68     1. Sets up logging
69     2. Checks for required dependencies
70     3. Creates and runs the main application
71     4. Handles any startup errors gracefully
72     """
73     try:
74         # Setup application logging
75         setup_logging()
76
77         # Check dependencies
78         if not check_dependencies():
79             logging.error("Dependency check failed. Exiting application.")
80             return
81
82         logging.info("All dependencies verified successfully")

```

```


83
84     # Create and run the main application
85     logging.info("Initializing main application...")
86     app = UnifiedVitalSignsApp()
87
88     logging.info("Starting application main loop...")
89     app.run()
90
91 except Exception as e:
92     error_msg = f"Fatal error during application startup: {e}"
93     logging.error(error_msg, exc_info=True)
94
95     # Show user-friendly error message
96     try:
97         messagebox.showerror(
98             "Application Error",
99             f"An error occurred while starting the application:\n\n{e}\n\nCheck the
100 log file for more details."
101         )
102     except:
103         # If even tkinter fails, print to console
104         print(f"FATAL ERROR: {error_msg}")
105
106     sys.exit(1)
107
108 finally:
109     logging.info("Vital Signs Monitor application shutdown complete")
110
111 if __name__ == "__main__":
112     """
113     Application entry point.
114
115     This ensures the application only runs when this file is executed directly,
116     not when imported as a module.
117     """
118     main()
119

```


Kode 33: Menjalankan GUI

5 Hasil

Sistem deteksi detak jantung dan pernapasan real-time telah berhasil diimplementasikan dengan antarmuka GUI yang menampilkan video feed langsung, kontrol monitoring, dan visualisasi sinyal secara simultan. Sistem dapat mendeteksi ROI pada area wajah untuk ekstraksi sinyal rPPG dan area bahu untuk deteksi pergerakan respirasi menggunakan MediaPipe. Dari hasil pengujian, sistem berhasil mengekstrak sinyal RGB dari ketiga channel warna dengan intensitas yang stabil dan konsisten. Setelah pemrosesan menggunakan metode POS dan filtering, sinyal rPPG yang dihasilkan menampilkan pola periodik yang jelas, menunjukkan bahwa sistem mampu mendeteksi variasi detak jantung dengan akurasi yang baik dalam rentang normal fisiologis.




(a) Afen-Hasil Sinyal rPPG




(b) Shintya-Hasil Sinyal rPPG

Gambar 2: Hasil Sinyal rPPG



(a) Afen-Hasil Sinyal Respirasi



(b) Shintya-Hasil Sinyal Respirasi

Gambar 3: Hasil Sinyal Respirasi

Untuk sinyal respirasi, sistem berhasil mendeteksi pergerakan vertikal bahu yang kemudian difilter menggunakan Butterworth band-pass filter untuk menghasilkan sinyal respirasi yang bersih dan stabil. Hasil pengukuran respiration rate menunjukkan konsistensi dalam rentang normal pernapasan manusia. Sistem dilengkapi dengan kontrol filter yang dapat disesuaikan dan mampu melakukan monitoring kontinyu dengan stabilitas sinyal yang baik. Meskipun sistem menunjukkan kinerja yang memuaskan untuk monitoring non-invasive menggunakan webcam standar, terdapat keterbatasan dalam hal sensitivitas terhadap kondisi pencahayaan dan memerlukan posisi subjek yang relatif stabil untuk hasil optimal, namun secara keseluruhan sistem dapat berfungsi sebagai alat monitoring vital sign yang efektif dan mudah digunakan.

6 Kesimpulan

Sistem deteksi detak jantung dan pernapasan real-time menggunakan webcam telah berhasil diimplementasikan dengan menggabungkan metode Plane Orthogonal-to-Skin (POS) untuk ekstraksi sinyal rPPG dan MediaPipe Pose untuk deteksi pergerakan respirasi. Sistem mampu melakukan monitoring simultan terhadap kedua parameter vital dengan antarmuka yang user-friendly, menampilkan visualisasi sinyal real-time, dan menghasilkan pengukuran yang berada dalam rentang normal fisiologis. Penerapan filtering Butterworth band-pass dan teknik denoising wavelet berhasil menghasilkan sinyal yang bersih dan stabil untuk analisis lebih lanjut. Meskipun terdapat keterbatasan terkait sensitivitas terhadap kondisi pencahayaan dan memerlukan posisi subjek yang stabil, sistem ini membuktikan potensi teknologi computer vision dan digital signal processing dalam aplikasi monitoring kesehatan non-invasive yang dapat diakses menggunakan perangkat sederhana seperti webcam, membuka peluang untuk pengembangan lebih lanjut dalam bidang telemedicine dan healthcare monitoring.

References

- [1] W. Wang, A. C. den Brinker, S. Stuijk, and G. de Haan, “Algorithmic principles of remote ppg,” *IEEE Transactions on Biomedical Engineering*, vol. 64, no. 7, pp. 1479–1491, 2017. [Online]. Available: <https://ieeexplore.ieee.org/document/7565547>
- [2] C. Massaroni, D. Lo Presti, D. Formica, S. Silvestri, and E. Schena, “Non-contact monitoring of breathing pattern and respiratory rate via rgb signal measurement,” *Sensors*, vol. 19, no. 12, 2019. [Online]. Available: <https://www.mdpi.com/1424-8220/19/12/2758>
- [3] SciPy Developers, “scipy.signal.butter — scipy v1.11.3 manual,” <https://docs.scipy.org/doc/scipy/reference/generated/scipy.signal.butter.html>, 2023, accessed: 2025-05-31.